



Med Software LLC

Created: 5/1/2011

Version: 1

Page: 1 of 1

Sperm Inventory

User Instruction



On the main toolbar in REsource, the **Sperm Inventory** option in the **Male** menu generates an inventory report for the patient selected. The inventory report contains any specimens currently in storage.

To generate and print an inventory report:

1 – Select the patient in the Patient Explorer and select **Sperm Inventory** from the **Male** menu

2 – The system will generate the report. Click the **Print** button in the upper left corner to print the document.

Inventory

Print...

100 %

1/1

Back

Forward

RMANJ

Reproductive Medicine Associates of New Jersey - Morristown
111 Madison Avenue, Suite 100
Morristown, NJ 07962

v1.4.1044.001

Tuesday, May 10, 2005

Inventory Summary

aaaa, bbbbbb T.

111-00-0173

Sperm Source	Frozen Date	Specimen ID	Origin	Check In	Tech In	Check Out	Outcome	Destination	Tech Out
b22222, a111.	12/15/2003	11410	Abington Memorial Hospital	12/15/2003		12/15/2003	n/a		
Guar Test 1, G1	12/01/2003	xxx	Abington Memorial Hospital	12/15/2003		12/15/2003	n/a		
b22222, a111.	12/15/2003	11412	Abington Memorial Hospital	12/15/2003		12/15/2003	n/a		
Guar Test 1, G1	12/15/2003	11413	Abington Memorial Hospital	12/15/2003		12/15/2003	n/a	Discard	
b22222, a111.	12/15/2003	11414	Abington Memorial Hospital	12/15/2003		12/15/2003	Transfer Out	Discard	
Wales, Joes M.	12/15/2003	11415	Abington Memorial Hospital	12/15/2003		02/18/2004	n/a	For Test	

TIP: There is also a **Sperm Inventory** report available in the **Reports** menu on the main toolbar. This report is an inventory of all sperm specimens stored in the cryo tanks and may take a few minutes to generate when clicked.