



3. Image boxes will come into view:



4. Move each image from the left-hand box to the right-hand box using the arrows in the middle. Highlight image to select it (blue arrow). This is moving the images into the Patient Images window (Patient menu - Images, or Ctrl G) .

Once they are all in there then you will go to the new stim and move them out of Pt Images:

1. Close the stim Images window
2. Go to the selected stim row you wish to move the images and create an ultrasound order, then open the Images window
3. Click the 'Move Img (ON)' button and move the ultrasound images from the right-hand box:



4. Repeat for each applicable stim date