


Embryo/Sperm Check In

- Added 'Store' button for Sperm/Embryo check in (red arrow). This button allows you to keep all the stored information from the previously stored sperm/oocyte/embryo check in.

Example:



Clicking 'Store' will keep 'Origin' and 'Specimen Type' info for the next specimen. However, a user will have to enter the new Specimen Id unless Pin Specimen Id box is checked off.

- Added Pin Specimen Id check box. If checked, the ID field will remain populated with value and will not clear after pressing 'Store'. (blue arrow)
- Cryo Tank Slots – when accessed from CheckIn
- Added Find Pat button. This will find first slot occupied by the selected patient's embryos/sperm.

Lab Order

The labs available list will now search incrementally as keys are typed.

Rx

Bug Fix

The med parameters screen was appearing when editing med type 'Apparatus' (i.e. Syringe) that has been added to a script but not yet saved.