


Patient Demo

On Resource Main Toolbar, click Patient, scroll to Demographics:



Country is now an editable droplist.

NOTE: Admin Setup: When reviewing the lkupCountries Table in the LookUp Table Manager: Fields/values entered that are not in the list, will appear in RED.

- The state droplist will repopulate with states when countries that have states (either USA, Mexico or Canada) is selected for country.
- Blank country defaults to USA states (like it always has).
- System switches were added to turn on hard validation for Country and State. When these switches are on, the demo form/patient details form will not save unless a valid State and country are entered or they are blank.
- Med Software have left these switches turned OFF and will turn them on at your request.

IVF Scheduler

Currently only 2 stim types (Cryo natural and synth) are seen in the system as Cryo start by the IVF Scheduler. A new attribute was added that can be configured to allow other stim types to be seen this way so they can be scheduled as Cryo starts in the scheduler. This might be useful if a new cryo stim type was added to the system.

LIS

Added new interpretive tolerance rule type. These rules may return a textual interpretation of a lab result based on a comparison of the result to a value or discreet ranges of values when the lab is posted. If an interpretation is returned, the lab comments field is set to this value and posted with the lab order result.

Example:

The screenshot shows a 'Lab Result' dialog box with the following fields and values:

- Lab Type: HIV 1/0/2
- Date Performed: 10/16/2011
- Performed by Third Party Lab: ☐
- Result Outcome: ☐ ? ☒ Normal ☐ Abnormal
- Result: <1.00
- FollowUp Date: (empty)
- Comment: NonReactive

A green arrow points from a text box to the 'Comment' field. The text box contains the following text:

This interpretation that accompanied the lab result '<1.00' was included in the comment field of this lab order

Buttons: OK, Cancel, Sign Off (checked)

Embryo/Sperm Check In

- Added 'Store' button for Sperm/Embryo check in (red arrow). This button allows you to keep all the stored information from the previously stored sperm/oocyte/embryo check in.

Example:

The screenshot shows a 'Check In' dialog box with the following fields and values:

- Specimen Source: (empty)
- Specimen ID: 12345
- Origin: Fresh
- Specimen Type: PatientEgg/PartnerSp
- Status: (empty)

Buttons: OK, Cancel, Store (highlighted with a red arrow), Pin Specimen Id (highlighted with a blue arrow).

Clicking 'Store' will keep 'Origin' and 'Specimen Type' info for the next specimen. However, a user will have to enter the new Specimen Id unless. Pin Specimen Id box is checked off.

- Added Pin Specimen Id check box. If checked, the ID field will remain populated with value and will not clear after pressing 'Store'. (blue arrow)
- Cryo Tank Slots – when accessed from CheckIn
- Added Find Pat button. This will find first slot occupied by the selected patient's embryos/sperm.

Rx

Bug Fix

The med parameters screen was appearing when editing med type 'Apparatus' (i.e. Syringe) that has been added to a script but not yet saved.